Double Flowmeter, DIN, O2, 0-15l/min



Manufacturer:	Greggersen
Product Number:	900745

Description

Double Flowmeter, O2, DIN

-PLUG-IN UNIT-

Medical double flowmeter,

pressure-compensated according to DIN EN ISO 15002,

with fine control valve for dosing the flow,

for connection to a terminal unit

Gas type: Oxygen / O2

Connection: Probe according to DIN 13260 Part 2

Outlet: 2 x 9/16-18 UNF

Flow rate: each 0...15 l/min continuously adjustable

Inlet pressure: 450 kPa +/- 50 kPa

Properties

© 2026, Greggersen Gasetechnik GmbH

Gas Type	O2-Oxygen

© 2026, Greggersen Gasetechnik GmbH