Air Motor Forano, Flush-Mounted/Cavity Wall



Manufacturer:	Greggersen
Product Number:	902090

Description

Terminal Unit Air Motor FORANO

For operating surgical tools with compressed air;
with exhaust air return; with non-return valve.
All-metal version - front plate and installation housing
made of stainless steel.

Installation type: Flush-mounted/cavity wall
Inlet: CU pipe 8x1 mm
Outlet: CU pipe 15x1 mm
Standard: DIN 13260-2 (Type 1), DIN EN ISO 9170-1

© 2026, Greggersen Gasetechnik GmbH