Compressed-Air Vacuum Regulator Specht -90 DIN Probe



Manufacturer:	Greggersen
Product Number:	900960

Description

Compressed-Air Vacuum Regulator Specht-90, Plug-in Unit DIN

Medical compressed-air vacuum regulator for vacuum generation,
with vacuum gauge, fine adjustment and quick-shut-off valve,
silencer and filter holder.
Technical Data:
Gas Type: Compressed Air / AIR
Inlet: Probe to DIN 13260 Part 2
Control Range: 0 to -85 kPa (0 to -0.85 bar)
Outlet: 9/16-18 UNF with vacuum hose barb

Properties

Gas Type	AIR Compressed Air

© 2026, Greggersen Gasetechnik GmbH