Vacuum Reg. Pirol 2.0 -30,VAC,PROBE,DIN



Manufacturer:	Greggersen
Product Number:	904606

Description

Vacuum Regulator Pirol -30, Plug-in Unit DIN

Medical vacuum regulator for vacuum dosing.
With 360° rotatable pressure gauge; large, ergonomic adjustment handwheel. The adjustment mechanism allows stepless, precise and intuitive control of the vacuum. The rotatable pressure gauge enables the user to read the set value from an optimal viewing angle.
The regulator unit meets the requirements of DIN EN ISO 10079-3 (high vacuum /high flow).
Material:
Housing: Aluminium, anodised
Handwheel: Plastic
Probe front part: Stainless steel
Technical data:

© 2026, Greggersen Gasetechnik GmbH

Gas type: Vacuum / VAC

Inlet: Probe according to DIN 13260 Part 2

Control range: 0 to -30 kPa ± 10%

Flow rate: > 20 l freeflow/min

Outlet: 9/16-18 UNF with vacuum tubing nipple

Properties

Gas Type	VAC-Vacuum

© 2026, Greggersen Gasetechnik GmbH