
Falke with Flow Paediatric O2 DIN Tubing Nipple



Manufacturer:	Greggersen
Product Number:	910202

Description

Falke with Flow, Paediatric,

Compact pressure regulator with , contents gauge;

incrementally adjustable flow rate.

Gas type: Oxygen / O2

Connection: Cylinder Connection for gas cylinders according to DIN 477, No. 9

Outlet: Tubing Nipple bottom

Flow Rate: 0...5 l/min.

Flow Settings: 0/ 0.1/ 0.3/ 0.5/ 0.6/ 0.8/ 1/ 2/ 3/ 4/ 5 l/min

Properties

Gas Type	O2-Oxygen

© 2026, Greggersen Gasetechnik GmbH