Airmotor Forano, Rail



Manufacturer:	Greggersen
Product Number:	902092

Description

Terminal Unit Air Motor FORANO

For operating surgical tools with compressed air;
with return line for exhaust air; with non-return valve.
All-metal version

Installation type: Rail
Inlet: CU pipe 8x1 mm
Outlet: CU pipe 15x1 mm
Standard: DIN 13260-2 (Type 1), DIN EN ISO 9170-1

Properties

Gas Type	AIR Compressed Air

© 2026, Greggersen Gasetechnik GmbH