Vacuum Regulator Spatz -16 DIN Probe



Manufacturer:	Greggersen
Product Number:	900899

Description

Vacuum Regulator Spatz-16, Plug-in Unit DIN

Medical vacuum regulator with vacuum gauge, fine adjustment and quick-acting shut-off valve.

Technical data:

Gas type: Vacuum / VAC

Inlet: Probe according to DIN 13260 Part 2

Control range: 0 to -16 kPa (0 to -160 mbar)

Outlet: 9/16-18 UNF with vacuum tubing nipple

© 2026, Greggersen Gasetechnik GmbH