Flowmeter, DIN, CO2, 0-12l/min



Manufacturer:	Greggersen
Product Number:	500011

Description

Flowmeter, Standard, CO2

Plug-in Unit

Medical flowmeter, pressure-compensated according to DIN EN ISO 15002, with fine adjustment valve for flow rate dosing, for coupling into a terminal unit

Technical Data:

Gas type: Carbon dioxide / CO2

Connection: Probe according to DIN 13260 Part 2

Outlet: 9/16-18 UNF

Flow rate: 0...12 l/min

Properties

Gas Type	CO2-Carbon Dioxide

© 2026, Greggersen Gasetechnik GmbH