# Vacuum Regulator High-Spatz -90 DIN Probe



| | |
|---|---|
| **Manufacturer:** | Greggersen |
| **Product Number:** | 900897 |

## Description

Vacuum Regulator High-Spatz-90, Plug-in Unit DIN

Medical vacuum regulator with vacuum gauge,
fine dosing and quick-action shut-off valve.
Technical data:
Gas type: Vacuum / VAC
Inlet: Probe to DIN 13260 Part 2
Control range: 0 to -85 kPa (0 to -0.85 bar)
Outlet: 9/16-18 UNF with vacuum tubing nipple

© 2026, Greggersen Gasetechnik GmbH