Habicht Rail, O2, NIST, Kolibri



Manufacturer:	Greggersen
Product Number:	902676

Description

Emergency Supply Unit

Suction and insufflation unit, for oxygen

RAIL-MOUNTED UNIT NIST

Combi-unit for attachment to an equipment rail; gas supply via a pressure regulator or a terminal unit; ejector regulator with vacuum gauge, fine control valve and silencer with filter holder; Kolibri compact flowmeter with 10 flow settings

Gas type: Oxygen / O2

Connection: NIST according to DIN EN ISO 18082

Flowmeter outlet: 9/16-18 UNF with tubing nipple

Ejector outlet: Vacuum tubing nipple

© 2026, Greggersen Gasetechnik GmbH

consisting of:

-Rail clamp with pressure inlet NIST

-Distribution block, spacer

-Kolibri standard 0...15 l/min

-Ejector regulator, 0 to -0.9 bar

© 2026, Greggersen Gasetechnik GmbH