
Falke with Standard Flow O2 DIN Tubing Nipple +1DIN



Manufacturer:	Greggersen
Product Number:	910204

Description

Falke with Flow, standard with 1 quick-release coupling,

Compact Pressure Regulator with contents gauge;

stepwise adjustable flow rate,

DIN quick-release coupling at a pressure of 450 kPa.

Gas type: Oxygen / O2

Connection: Cylinder Connection for gas cylinders according to DIN 477, No. 9

Outlet: Tubing Nipple / 1 quick-release coupling according to DIN 13260 Part 2

Output: 0...15 l/min Tubing Nipple - 120 l/min quick-release coupling

Flow Settings: 0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Properties

Gas Type	O2-Oxygen

© 2026, Greggersen Gasetechnik GmbH