Nebuliser Plug-in Valve with Clamp O2



Manufacturer:	Greggersen
Product Number:	900776

Description

Fitis Nebuliser Plug-in Valve with Clamp O2

Plug-in Unit

Nebuliser plug-in valve, preset with side clamp;

for coupling into a terminal unit

Technical Data:

Gas type: Oxygen

Connection: Probe to DIN 13260 Part 2

Outlet: 9/16-18 UNF/Tubing Nipple

Flow Rate: 5 l/min ± 20%

Properties

Gas Type	O2-Oxygen

© 2026, Greggersen Gasetechnik GmbH