Compressed-Air Vacuum Regulator Specht -10 Rail NIST



Manufacturer:	Greggersen
Product Number:	902672

Description

Compressed-Air Vacuum Regulator Specht-10, Rail-Mounted Unit NIST

Medical compressed-air vacuum regulator for vacuum generation,
with vacuum gauge, fine adjustment and quick-acting valve,
silencer and filter holder.
For mounting to an equipment rail (25x10 mm).
The gas inlet is integrated in the rail clamp (NIST).

Accessories required for operation:
Low-Pressure Hose AIR with NIST outlet
(not included in the scope of delivery)
Technical data:
Gas type: Compressed air / AIR

© 2026, Greggersen Gasetechnik GmbH

Inlet: NIST housing according to DIN EN ISO 18082

Control range: 0 to -10 kPa (0 to -100 mbar)

Outlet: 9/16-18 UNF with vacuum tubing nipple

Properties

Gas Type	AIR Compressed Air

© 2026, Greggersen Gasetechnik GmbH