
Forano, Front Part, DIN, O2, ISO



Manufacturer:	Greggersen
Product Number:	905100

Description

Terminal Unit FORANO - Front Part DIN

ISO 32

Gas type: Oxygen / O2

Inscription: oxygen / O2

Gas-specific socket assembly for installation

on base part FORANO.

Socket assembly in accordance with DIN 13260 Part 2

with parked and operating position. Unlocking via release ring.

Colour coding:

according to DIN EN ISO 9170-1 (ISO32 - white)

Release ring inscription: printed

Language: English

Cartridge incl. sealing elements

© 2026, Greggersen Gasetechnik GmbH

Technical data:

Gas type: Oxygen / O2

Inscription: oxygen / O2

Colour coding: according to ISO 9170-1 (white)

Front diameter: Ø 42mm

Scope of delivery:

- gas-specific socket assembly
- release ring incl. printed inscription
- cartridge with sealing elements

Properties

Gas Type	O2-Oxygen

© 2026, Greggersen Gasetechnik GmbH