Compressed-Air Vacuum Regulator Specht -90 Rail Anti



Manufacturer:	Greggersen
Product Number:	500249

Description

Compressed-Air Vacuum Regulator Specht-90, Rail-Mounted Unit NIST, non-magnetic

Medical compressed-air vacuum regulator for vacuum generation,
with vacuum gauge, fine adjustment and quick-acting shut-off valve,
silencer and filter holder.
For attachment to an equipment rail (25x10 mm).
The gas inlet is integrated into the rail clamp (NIST).

Accessories required for operation:
Low-Pressure Hose AIR with NIST outlet
(not included in the scope of delivery)
Technical data:
Gas type: Compressed air / AIR

© 2026, Greggersen Gasetechnik GmbH

Inlet: NIST housing according to DIN EN ISO 18082

Control range: 0 to -85 kPa (0 to -0.85 bar)

Outlet: 9/16-18 UNF with vacuum tubing nipple

Properties

Gas Type	AIR Compressed Air

© 2026, Greggersen Gasetechnik GmbH