Habicht Rail, O2, NIST, Kolibri+Coupling



Manufacturer:	Greggersen
Product Number:	902677

Description

Emergency Supply Unit

Suction and insufflation unit, for oxygen

RAIL-MOUNTED UNIT NIST

Combi-unit for attachment to an equipment rail

Gas supply via a pressure regulator or

a gas terminal unit; ejector regulator with

vacuum gauge, fine control valve and silencer

with filter mount, Kolibri compact flowmeter with

10 flow settings and additional 5 bar outlet

as central gas supply coupling.

Gas type: Oxygen / O2

Connection: NIST according to DIN EN ISO 18082

consisting of:

© 2026, Greggersen Gasetechnik GmbH

-Rail clamp with NIST pressure inlet

-Distribution block, spacer

-Kolibri Standard 0...15 l/min

-Ejector regulator 0 to -0.9 bar

-Quick-connect coupling according to DIN 13260 Part 2

© 2026, Greggersen Gasetechnik GmbH