
Falke with FLOW Standard O2 DIN Tubing Nipple



Manufacturer:	Greggersen
Product Number:	910200

Description

Falke with Flow, Standard,

Compact pressure regulator with contents gauge;

step-adjustable flow rate.

Gas type: Oxygen / O2

Connection: Cylinder connection for gas cylinders according to DIN 477, no. 9

Outlet: Tubing nipple below

Flow rate: 0...15 l/min

Flow settings: 0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 l/min

Properties

Gas Type	O2-Oxygen

© 2026, Greggersen Gasetechnik GmbH