
Angled Probe CO2 Hose Barb 7mm, DIN,ISO V4A



Manufacturer:	Greggersen
Product Number:	505154

Description

Angled Probe for Carbon Dioxide

with 7 mm Hose Barb and Non-Return Valve

Probe Front Part made of 1.4404 (Stainless Steel)

Gas Type: Carbon Dioxide / CO2

Connection: Probe to DIN 13260 Part 2

Outlet: 7 mm Hose Barb

Colour Coding: to ISO 9170-1 (grey)

Properties

Gas Type	CO2-Carbon Dioxide

© 2026, Greggersen Gasetechnik GmbH