Vacuum Reg. Pirol 2.0 -30,VAC,RAIL,NIST



Manufacturer:	Greggersen
Product Number:	904607

Description

Vacuum Regulator Pirol -30, Rail-Mounted Unit NIST

Medical vacuum regulator for dosing vacuum.
With 360° rotatable gauge; large, ergonomic
adjustment handwheel. The adjustment mechanism enables continuous,
precise and intuitive vacuum regulation. The rotatable
gauge allows the user to read the
set value from an optimal viewing angle.
The regulator unit meets the requirements of
DIN EN ISO 10079-3 (high vacuum / high flow).
Material:
Housing: aluminium, anodised
Handwheel: plastic
Rail clamp: aluminium, anodised
Technical data:

© 2026, Greggersen Gasetechnik GmbH

Gas type: Vacuum / VAC

Inlet: NIST housing according to DIN EN ISO 18082

Control range: 0 to -30 kPa ± 10%

Flow rate: > 20 l freeflow/min

Outlet: 9/16-18 UNF with vacuum tubing nipple

Properties

Gas Type	VAC-Vacuum

© 2026, Greggersen Gasetechnik GmbH