Double Flowmeter, DIN, AIR, 0-15l/min



Manufacturer:	Greggersen
Product Number:	900748

Description

Double Flowmeter

-PLUG-IN UNIT-

Medical double flowmeter,

pressure-compensated according to DIN EN ISO 15002,

with fine control valve for flow rate adjustment,

for connection to a terminal unit

Gas type: Compressed Air

Connection: Probe according to DIN 13260 Part 2

Outlet: 2 x 9/16-18 UNF

Flow rate: each 0 ... 15 l/min continuously adjustable

Inlet pressure: 450 kPa +/- 50 kPa

Properties

© 2026, Greggersen Gasetechnik GmbH

Gas Type	AIR Compressed Air

© 2026, Greggersen Gasetechnik GmbH