Portable Emergency Supply Habicht Tube



Manufacturer:	Greggersen
Product Number:	902696

Description

Emergency supply unit,
portable, without oxygen cylinder
- with tube flowmeter -

Emergency supply unit, oxygen flowmeter
and ejector for vacuum generation
consisting of:
Stainless steel carrying frame with stand, to accommodate
a 3 l O2 cylinder, suitable
for mounting on equipment rail or bed,
central supply block with NIST inlet,
consisting of:
Flowmeter with fine dosing valve for regulating the flow
and ejector with vacuum gauge and fine dosing valve;

Gas type: Oxygen / O2
Connection: NIST in accordance with DIN EN ISO 18082
Outlet: Tubing nipple / vacuum tubing nipple
Performance: FLOW : 0...15 l/min

© 2026, Greggersen Gasetechnik GmbH

VAC: adjustable between

0 and -85 kPa (0 to -0.85 bar)

© 2026, Greggersen Gasetechnik GmbH