# Vacuum Regulator Spatz -16 DIN Probe long




| | |
|---|---|
| **Manufacturer:** | Greggersen |
| **Product Number:** | 900898 |

## Description

Vacuum Regulator Spatz-16, Plug-in Unit DIN extended version (+57 mm)

Medical vacuum regulator with vacuum gauge, fine dosing and quick-shut-off valve.
Technical data:
Gas type: Vacuum / VAC
Inlet: Probe according to DIN 13260 Part 2
Control range: 0 to -16 kPa (0 to -160 mbar)
Outlet: 9/16-18 UNF with vacuum tubing nipple

© 2026, Greggersen Gasetechnik GmbH