Vacuum Regulator Specht -16 AIR DIN Probe



Manufacturer:	Greggersen
Product Number:	900958

Description

Compressed-Air Vacuum Regulator Specht-16, Plug-in Unit DIN

Medical compressed-air vacuum regulator for vacuum generation,
with vacuum gauge, fine dosing and quick-action shut-off valve,
silencer and filter mount.
Technical data:
Gas type: Compressed air / AIR
Inlet: Probe according to DIN 13260 Part 2
Control range: 0 to -16 kPa (0 to -160 mbar)
Outlet: 9/16-18 UNF with vacuum tubing nipple

Properties

Gas Type	AIR Compressed Air

© 2026, Greggersen Gasetechnik GmbH