
Forano, Front Part, DIN, AIR ,ISO



Manufacturer:	Greggersen
Product Number:	905101

Description

Terminal Unit FORANO - Front Part DIN

ISO 32

Gas type: Compressed Air / AIR

Labelling: medical air / AIR

Gas-specific socket assembly for mounting

on FORANO base part.

Socket assembly according to DIN 13260 Part 2

with parking and operating position.

Release via release ring.

Colour coding:

according to DIN EN ISO 9170-1 (ISO32 - black/white)

Release ring labelling: printed

Language: English

© 2026, Greggersen Gasetechnik GmbH

Cartridge incl. sealing elements

Technical data:

Gas type: Compressed Air / AIR

Labelling: medical air / AIR

Colour coding: according to ISO 9170-1 (black/white)

Front diameter: Ø 42mm

Scope of delivery:

- gas-specific socket assembly
- release ring incl. printed labelling
- cartridge with sealing elements

Properties

Gas Type	AIR Compressed Air

© 2026, Greggersen Gasetechnik GmbH