
Ventus 22 mm Valve Insert AIR-800



Manufacturer:	Greggersen
Product Number:	900787

Description

Valve insert for Ventus 22 mm

Gas type : AIR-800

Design complies with DIN EN ISO 7396-1

Valve assembly consisting of:

- Shut-off valve DN 20, G3/4"
- CU pipe 22mm inlet / outlet
- Pressure gauge 0...16 bar
- Pressure sensor 0...16
- NIST emergency feed point
- Physical separation point
- Gas type labelling

© 2026, Greggersen Gasetechnik GmbH

Properties

Fulfilled Standard	DIN EN ISO 7396-1
Gas Type	AIR-800-Compressed Air 8bar
NIST Emergency Inlet Point	1 NIST Emergency Inlet Point

© 2026, Greggersen Gasetechnik GmbH